

Purpose: To create simple galvanic cells to measure and compare actual cell potentials to theoretical cell potentials.

Procedures:

PART ONE

$\text{Zn(s)} / \text{Zn}^{2+}(0.05\text{M})$

$\text{Cu}^{2+}(0.05\text{M}) / \text{Cu(s)}$

$\text{Zn(s)} + \text{Cu}^{2+}_{(\text{aq})} \rightarrow \text{Zn}^{2+}_{(\text{aq})} + \text{Cu(s)}$

- 1) clean a strip of each metal, Cu(s) and Zn(s) with emery paper.
- 2) clean and rinse two small beakers with water.
- 3) Fill one beaker with 20 mL of copper sulfate and the other with 20 mL of zinc sulfate.
- 4) connect the two beakers with a salt bridge which is a 1.0 M KNO_3 solution in water.
- 5) Place the copper in the copper sulfate solution and the zinc metal in the zinc sulfate solution.
- 6) complete the circuit by connecting the voltmeter ~~to the~~.

7) Record the reading on the voltmeter.

CELL DIAGRAM

$\text{Zn(s)} | \text{Zn}^{2+}(\text{aq}, 0.05\text{M}) || \text{Cu}^{2+}(\text{aq}, 0.05\text{M}) | \text{Cu(s)}$

PART TWO

- Refer to question # 2 in Prelab.

CELL DIAGRAM

$\text{Cu}^{2+}(\text{aq}, 0.05\text{M}) | \text{Cu(s)} || \text{Cu}^{2+}(0.00005\text{M}) | \text{Cu(s)}$

PART THREE

- ① obtain an E_{cell} value from the TA and compare this value

to the standard potential values in prelab question # 1.

- ② use the cell potential that is the closest higher value than the value obtained from the TA. ~~since our value was 0.97, we found that~~
- ③ since the electrochemical cell including the oxidizing agent Cu^{2+} and reducing agent Zn^{2+} has a value of 1.03 V, which is closest to and greater than the obtained value of 0.97, use this value.
- ④ Use the Nernst equation to solve for Q , using E_{cell} and the standard potential as follows:

$$0.97\text{V} = 1.103\text{V} - (0.0592/2) \log Q$$

$$E_{\text{cell}} = E^{\circ}_{\text{cell}} - \frac{0.0592}{2} \log Q$$

$$\log Q = \frac{2 \times (0.97 - 1.103)}{-0.0592}$$

$$= 4.4932$$

$$Q = 311345.9661$$

$$= 0.059$$

- ⑤ since the highest value of Zn^{2+} is 0.05 M, and $Q = [\text{Zn}^{2+}_{\text{aq}}] / [\text{Cu}^{2+}_{\text{aq}}]$, solve for the concentration of Q .

$$311345.9661 = \frac{0.05\text{M}}{[\text{Cu}^{2+}]}$$

$$[\text{Cu}^{2+}] = 1.6 \times 10^{-6} \text{ M}$$

SIGNATURE

DATE

WITNESS/TA

DATE