Frelab Questions

Melling Point of Asprin: 135°C Source Haynes, W , Lide, D. Cic Handbeck of Chemistry and Physics, 90 in ed Florida: (rc press, 2009 (3-8)

## Furpose:

This purpose of this experiment 12 to perform several experiments utilizing the asprin product oblained from experiment # 4.

## Ficcodurs

The experiment was carried out as described in the 2011/2012 Chem 1403/14A3 Lab Manual, organic Chemistry - Z.

## CISSEKVATIONS


Grant gilld calculation 1. YILIN = Actual x 100 Theoretical

$$= \frac{0.8791}{2.609} \times 100 = 33.77.$$

6 see discussion for explanation.

Melling Point: 124-1260C. oses discussion por explanation.

## Z. Thin Layer chiemategraph (ILC) TA Demonstration



SA = = alget salreyelle acid (starting material) r = product ASA ASA = Aspirin.

6 sec Discussion for explanation

Reactions of Aspirin 4. b) Cossivations (aughtative)

Test rubs 1: Aspirin and Natl ( 3 observation: Eubbles formed on top of solution and almost all aspirin crystals dissolved. Leaves a clear, water-like solution.

Tost tubs 2: salicyclic Acid and Nations observation Thick layer of bubbles formed immediatily followed by the cornation of a white percipitati which ramains d in solution.

Test Tube 3: Benzyl Alcohol and NAHCOS Observation: No chemical reaction occured. Remained separate from test tube one solution floated on top of the other.

Gase discussion for analysis of in dividual reactions.

SIGNATURE	DATE	WITNESS/TA	DATE
			TOWNER LINDER CORV CHEET REFORE WRITING