LAB PARTNER

Discussion: Nach (ag) + Heliag) - Nacliag) + Hz (12) During this experiment, an Nach solution was titrated against an HCI solution. NaOH was released into the MCI solution until it had completely reacted with the HCI, er basically when the equivalence point of the reaction was reached. This point in the reaction was Visible as a faint pink coleured solution, as a result of the phenolothelein indicator which was pieced in the Hel sciution, prior to the start of the titration. At this point in the reaction. the moles of nci and Nach are equal to each other, since the balanced enemical reaction, the volumes of NacH, and the standardized concentration of tici are known, the concentration of Ngoti was egiculated To obtain a more scientifically sound results, the experiment was conducted three times, and the mean of Nach concentration was calculated and

it is also important to take into consideration the sources of error that would affect results. Specifically, it was difficult to obtain a precise volume of exactly 10.00 m L of HCl solution. Furthermore, due to the nature of phenolphalein which turns pink when the solution is slighty bask,

exact equivalence point of the thration, and there was a slight excess of Nach in the solution.

11

recorded.